

Product Datasheet

CD10 antibody (APC) (orb432615)

Description	Mouse monoclonal antibody against CD10 conjugated to APC
Species/Host	Mouse
Reactivity	Human
Conjugation	APC
Tested Applications	FC
Immunogen	Partially purified cell membrane antigens from fresh leukemia cells
Target	CD10
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin5% sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to allophycocyanin (APC) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.