

Product Datasheet

CD71 antibody (orb432633)

Description	Mouse monoclonal antibody to CD71
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	PHA activated rat lymphocytes.
Target	CD71
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.