

Product Datasheet

TNF ALPHA antibody (RPE) (orb432638)

Description	Rat monoclonal antibody against TNF ALPHA conjugated to RPE
Species/Host	Rat
Reactivity	Human
Conjugation	RPE
Tested Applications	FC
Target	TNF ALPHA
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin5% sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.