

Product Datasheet

MHC CLASS I RT1Ac antibody (orb432654)

Description	Mouse monoclonal antibody to MHC CLASS I RT1Ac
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr
Immunogen	PHA activated rat lymphocytes.
Target	MHC CLASS I RT1Ac
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.