

Product Datasheet

CD40 antibody (orb432674)

Description	Mouse monoclonal antibody to CD40
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IHC-P, IP
Immunogen	CD40 Ig(Fc) fusion protein containing the EC region of human CD40 and Fc region of human IgG.
Target	CD40
Preservatives	0.1% sodium azide (NaN ₃). Tris buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.