

Product Datasheet

CD45RB antibody (FITC) (orb432707)

Description	Mouse monoclonal antibody against CD45RB conjugated to FITC
Species/Host	Mouse
Reactivity	Bovine
Conjugation	FITC
Tested Applications	FC
Immunogen	Bovine thymocytes
Target	CD45RB
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.