

Product Datasheet

BRCA1 antibody (orb432724)

Description	Mouse monoclonal antibody to BRCA1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, IP, WB
Immunogen	Recombinant protein corresponding to the N-Terminal region of human BRCA1
Target	BRCA1
Preservatives	0.1% sodium azide (NaN ₃). TRIS-glycine buffered saline, NaCl
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(N-TERMINAL)
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.