

Product Datasheet

CD8 BETA antibody (orb432744)

Description	Mouse monoclonal antibody to CD8 BETA
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP, WB
Immunogen	Normal human T lymphocytes.
Target	CD8 BETA
Preservatives	1% bovine serum albumin. Phosphate buffered saline
Form/Appearance	Purified IgG - lyophilized
Concentration	1.0 mg/ml
Storage	-20°C after reconstitution
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.