

Product Datasheet

CD44v4 antibody (FITC) (orb432750)

Description	Mouse monoclonal antibody against CD44v4 conjugated to FITC
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Immunogen	Glutathione S Transferase (GST) fusion protein corresponding to the variable domains (v3 to v10) of human CD44.
Target	CD44v4
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin. Tris buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.