

Product Datasheet

CD31 antibody (orb432786)

Description	Mouse monoclonal antibody to CD31
Species/Host	Mouse
Reactivity	Porcine
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr
Immunogen	Porcine CD31/human IgGFc fusion protein.
Target	CD31
Preservatives	0.09% sodium azide (NaN ₃)
Form/Appearance	Tissue culture supernatant - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
Purity	S/N
Expiration Date	12 months from date of receipt.