

Product Datasheet

CD4 ALPHA antibody (FITC) (orb432788)

Description	Mouse monoclonal antibody against CD4 ALPHA conjugated to FITC
Species/Host	Mouse
Reactivity	Porcine
Conjugation	FITC
Tested Applications	FC
Immunogen	Leucocytes isolated from porcine gut lamina propria.
Target	CD4 ALPHA
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.