

Product Datasheet

CD8 antibody (orb432814)

Description	Rat monoclonal antibody to CD8
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr
Target	CD8
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Preservative Free
Expiration Date	12 months from date of receipt.