

Product Datasheet

CD11c antibody (orb432821)

Description	Mouse monoclonal antibody to CD11c
Species/Host	Mouse
Reactivity	Canine
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IP
Target	CD11c
Preservatives	0.1% sodium azide
Form/Appearance	Tissue culture supernatant - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	S/N
Expiration Date	12 months from date of receipt.