

Product Datasheet

IFN GAMMA antibody (Biotin) (orb432828)

Description	Mouse monoclonal antibody against IFN GAMMA conjugated to Biotin
Species/Host	Mouse
Reactivity	Bovine
Conjugation	Biotin
Tested Applications	ELISA
Target	IFN GAMMA
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to biotin - liquid
Concentration	0.5 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.