

Product Datasheet

EGF R antibody (orb432831)

Description	Rat monoclonal antibody to EGF R
Species/Host	Rat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, IHC-Fr, IHC-P, IP
Immunogen	Extracellular domain of human EGF-receptor from head and neck carcinoma.
Target	EGF R
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.