

Product Datasheet

ESTROGEN RECEPTOR ALPHA antibody (orb432837)

Description	Mouse monoclonal antibody to ESTROGEN RECEPTOR ALPHA
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IH-Fr, IHC-P, WB
Immunogen	Recombinant human estrogen receptor (alpha form).
Target	ESTROGEN RECEPTOR ALPHA
Preservatives	0.09% sodium azide (NaN ₃)
Form/Appearance	Tissue culture supernatant - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	S/N
Expiration Date	12 months from date of receipt.