

Product Datasheet

CD21 antibody (orb432844)

Description	Mouse monoclonal antibody to CD21
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IH-Fr, IHC-P
Immunogen	Recombinant human CD21.
Target	CD21
Preservatives	0.09% sodium azide (NaN ₃)
Form/Appearance	Tissue culture supernatant - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	S/N
Expiration Date	12 months from date of receipt.