

Product Datasheet

CD8 antibody (orb432852)

Description	Mouse monoclonal antibody to CD8
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IH-Fr, IHC-P, WB
Immunogen	Synthetic peptide derived from the carboxy terminal region of the human CD8 alpha chain coupled to a N-terminal cysteine, with the sequence C-KSDGKPSLSARYV. The peptide was coupled to bovine serum albumin and keyhole limpet hemocyanin.
Target	CD8
Preservatives	0.09% sodium azide (NaN3)
Form/Appearance	Tissue culture supernatant - liquid
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	S/N
Expiration Date	12 months from date of receipt.