

Product Datasheet

CD81 antibody (FITC) (orb432892)

Description	Mouse monoclonal antibody against CD81 conjugated to FITC
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Immunogen	OCI-LY8 cells aggregated by 5A6 (another CD81 antibody)
Target	CD81
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.