

Product Datasheet

CD81 antibody (orb432895)

Description	Mouse monoclonal antibody to CD81
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, IHC-P, IP, WB
Immunogen	OCI-LY8 cells aggregated by 5A6 (another CD81 antibody)
Target	CD81
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.