

Product Datasheet

MARCO antibody (orb432904)

Description	Rat monoclonal antibody to MARCO
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, IHC-Fr, WB
Immunogen	CHO cells expressing mouse MARCO.
Target	MARCO
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.