

Product Datasheet

CD163 antibody (orb432914)

Description	Mouse monoclonal antibody to CD163
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	Leucocytes harvested from the pleural cavity of patients with idiopathic spontaneous pneumothorax
Target	CD163
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.