

Product Datasheet

M13 BACTERIOPHAGE antibody (Biotin) (orb432929)

Description	Mouse monoclonal antibody against M13 BACTERIOPHAGE conjugated to Biotin
Species/Host	Mouse
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA, FC
Immunogen	M13 Phage coat proteins.
Target	M13 BACTERIOPHAGE
Preservatives	0.08% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to biotin - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG3
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.