

Product Datasheet

CD147 antibody (APC) (orb432948)

Description	Mouse monoclonal antibody against CD147 conjugated to APC
Species/Host	Mouse
Reactivity	Human
Conjugation	APC
Tested Applications	FC
Immunogen	CD147Rg, consisting of the cDNA coding for the extracellular region of human CD147.
Target	CD147
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin5% sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to allophycocyanin (APC) - lyophilized
Concentration	0.1mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.