

Product Datasheet

CD252 antibody (orb432960)

Description	Rat monoclonal antibody to CD252
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr
Immunogen	Rat CD4 and mouse CD134 ligand Chimaeric protein.
Target	CD252
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.