

Product Datasheet

FELINE LEUKAEMIA VIRUS gp70 antibody (orb432965)

Description	Mouse monoclonal antibody to FELINE LEUKAEMIA VIRUS gp70
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IHC-P, IP, WB
Target	FELINE LEUKAEMIA VIRUS gp70
Preservatives	0.1% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.