

Product Datasheet

hCG antibody (orb432975)

Description	Mouse monoclonal antibody to hCG
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IA, IHC-P, RIA
Immunogen	hCG.
Target	hCG
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Application notes	(BETA 1 EPITOPE)
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.