

Product Datasheet

GM-CSF antibody (orb432983)

Description	Mouse monoclonal antibody to GM-CSF
Species/Host	Mouse
Reactivity	Sheep
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	CHO expressed recombinant ovine GM-CSF.
Target	GM-CSF
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.