

Product Datasheet

CD23 antibody (FITC) (orb432992)

Description	Mouse monoclonal antibody against CD23 conjugated to FITC
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Immunogen	Immunoprecipitated CD23 antigen from purified plasma membranes from EBV-transformed human lymphoblastoid B cell lines
Target	CD23
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.