


Product Datasheet

KDM5C antibody (orb43302)

Description	KDM5C antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Lysine-specific demethylase 5C protein (45-324AA)
Target	KDM5C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P41229
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MCF7 cell...