

Product Datasheet

CD137 antibody (orb433020)

Description	Rat monoclonal antibody to CD137
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Murine CD137 human Fc fusion protein
Target	CD137
Preservatives	0.1% sodium azide (NaN ₃). TRIS buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.