

Product Datasheet

CD68 antibody (orb433031)

Description	Rat monoclonal antibody to CD68
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, IHC-Fr, IHC-P, IP, WB
Immunogen	Purified Concanavalin A acceptor glycoprotein from P815 cell line.
Target	CD68
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline.
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.