

Product Datasheet

CD200 antibody (Biotin) (orb433033)

Description	Rat monoclonal antibody against CD200 conjugated to Biotin
Species/Host	Rat
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	FC
Immunogen	Mouse CD200-rat CD4 fusion protein.
Target	CD200
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to biotin - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.