

Product Datasheet

CD200 RECEPTOR 1 antibody (RPE) (orb433046)

Description	Mouse monoclonal antibody against CD200 RECEPTOR 1 conjugated to RPE
Species/Host	Mouse
Reactivity	Rat
Conjugation	RPE
Tested Applications	FC
Immunogen	Membrane fraction of thioglycollate-elicited rat peripheral cells.
Target	CD200 RECEPTOR 1
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin5% sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.