

Product Datasheet

CD18a antibody (RPE) (orb433073)

Description	Mouse monoclonal antibody against CD18a conjugated to RPE
Species/Host	Mouse
Reactivity	Porcine
Conjugation	RPE
Tested Applications	FC
Immunogen	Porcine large granular lymphocytes.
Target	CD18a
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin5% sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.