

Product Datasheet

CD240 antibody (orb433079)

Description	Mouse monoclonal antibody to CD240
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP
Immunogen	Human erythrocytes.
Target	CD240
Preservatives	0.1% sodium azide (NaN ₃). TRIS buffered glycine
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.