

Product Datasheet

CD236 antibody (orb433085)

Description	Mouse monoclonal antibody to CD236
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, Hemagglutination, IH-Fr, WB
Immunogen	Human erythrocytes.
Target	CD236
Preservatives	0.1% sodium azide (NaN ₃) ≤ 100mM glycine. TRIS buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.5mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.