

Product Datasheet

CD49b antibody (orb433100)

Description	Mouse monoclonal antibody to CD49b
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IP
Immunogen	Purified human b1 integrin preparation from HT1080 fibrosarcoma cell extract.
Target	CD49b
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.