

Product Datasheet

CD49f antibody (orb433106)

Description	Mouse monoclonal antibody to CD49f
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, IHC-Fr
Immunogen	Rat oesophageal carcinoma cell line B2T.
Target	CD49f
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.