

Product Datasheet

CD172a antibody (orb433114)

Description	Mouse monoclonal antibody to CD172a
Species/Host	Mouse
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Bovine afferent veiled dendritic cells
Target	CD172a
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.