

Product Datasheet

TUBULIN BETA 3 antibody (FITC) (orb433128)

Description	Mouse monoclonal antibody against TUBULIN BETA 3 conjugated to FITC
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	Synthetic peptide, ESESQGPK, corresponding to amino acids 441-448 of human class III beta tubulin coupled to Keyhole Limpet Hemocyanin (KLH). This sequence is widely conserved across species.
Target	TUBULIN BETA 3
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.