

Product Datasheet

CD1w2 antibody (orb433139)

Description	Mouse monoclonal antibody to CD1w2
Species/Host	Mouse
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IP
Target	CD1w2
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.