

Product Datasheet

CD1w2 antibody (RPE) (orb433140)

Description	Mouse monoclonal antibody against CD1w2 conjugated to RPE
Species/Host	Mouse
Reactivity	Bovine
Conjugation	RPE
Tested Applications	FC
Target	CD1w2
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin5% sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.