

Product Datasheet

CD284 antibody (Biotin) (orb433147)

Description	Mouse monoclonal antibody against CD284 conjugated to Biotin
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	FC
Immunogen	Ba/F3 cell line expressing TLR4 (CD284).
Target	CD284
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to biotin - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.