

Product Datasheet

CD284 antibody (orb433148)

Description	Mouse monoclonal antibody to CD284
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FACS, IP, WB
Immunogen	Ba/F3 cell line expressing TLR4 (CD284).
Target	CD284
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Low Endotoxin
Expiration Date	12 months from date of receipt.