

## Product Datasheet

### CD284 antibody (FITC) (orb433150)

|                            |                                                                                          |
|----------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody against CD284 conjugated to FITC                               |
| <b>Species/Host</b>        | Mouse                                                                                    |
| <b>Reactivity</b>          | Human                                                                                    |
| <b>Conjugation</b>         | FITC                                                                                     |
| <b>Tested Applications</b> | FC                                                                                       |
| <b>Immunogen</b>           | Ba/F3 cell line expressing TLR4 (CD284).                                                 |
| <b>Target</b>              | CD284                                                                                    |
| <b>Preservatives</b>       | 0.09% sodium azide (NaN <sub>3</sub> )1% bovine serum albumin. Phosphate buffered saline |
| <b>Form/Appearance</b>     | Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid           |
| <b>Concentration</b>       | 0.1 mg/ml                                                                                |
| <b>Storage</b>             | +4°C, -20°C if preferred                                                                 |
| <b>Note</b>                | For research use only                                                                    |
| <b>Isotype</b>             | IgG2a                                                                                    |
| <b>Clonality</b>           | Monoclonal                                                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                          |