

Product Datasheet

CD120a antibody (RPE) (orb433165)

Description	Rat monoclonal antibody against CD120a conjugated to RPE
Species/Host	Rat
Reactivity	Mouse
Conjugation	RPE
Tested Applications	FC
Target	CD120a
Preservatives	0.1% sodium azide (NaN ₃)0.5% bovine serum albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - liquid
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.