

Product Datasheet

CD18 antibody (orb433170)

Description	Mouse monoclonal antibody to CD18
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP, WB
Immunogen	Peripheral blood mononuclear cells.
Target	CD18
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(ACTIVATION EPITOPE)
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.