

Product Datasheet

CD18 antibody (RPE) (orb433172)

Description	Mouse monoclonal antibody against CD18 conjugated to RPE
Species/Host	Mouse
Reactivity	Human
Conjugation	RPE
Tested Applications	FC
Immunogen	Peripheral blood mononuclear cells.
Target	CD18
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Application notes	(ACTIVATION EPITOPE)
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.