

Product Datasheet

HLA A2 antibody (orb433186)

Description	Mouse anti-Human HLA A2 antibody, clone BB7.2 recognizes the human HLA-A2 histocompatibility antigen...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FACS, IHC-Fr, IP
Immunogen	Papain solubilized HLA-A2
Target	HLA A2
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Low Endotoxin
Expiration Date	12 months from date of receipt.